

# PROCESSOR QUALITY CONTROL

Sensitometer and Densitometer solution from AFD Technologies, Inc.

## Processor Quality Control is still required.

They're never gone away, but recently there have been an increase in inspections by many state agencies. However, as we move quickly towards an all digital world, it is harder to find quality sensitometers and densitometers. Many manufacturers have simply stopped making these products. Used equipment is, quite frankly, pretty beat up.

## We have the QA products you need.

We have never stopped providing good QA products, including sensitometers and densitometers. In fact, we have added new models to our offerings.

Shown below is our most popular housing. It's a combination sensitometer and densitometer. With our Densonorm, you always have both halves of the "operation" at your fingertips. It's one piece of equipment to keep track of, not two.

Here are the three most popular models (if you are a European customer, you may need a tracable calibrated unit, in which case call us).

Densonorm 21-ECO - a combination of a 21-step sensitometer and a point densitometer.

Densonorm 21-E - automatically calculates Relative Speed and Average Gradient after you have manually measured each step.

Densonorm 21i - adds an automatic film reader to the "E" version and also calculates all of the important measurements.

Each unit produces the information needed for daily quality control, but with the "21i" all of the work is done seconds after the unit reads the film strip. This saves time and eliminates the hassle of film QA.

**Step 2 "I":** With the "21i" slip film strip in here and the unit will pull it through its reader and make all calculations.

**Step 1:** Produce sensitometry strip here.



**Step 4:** Read out the results here.

**Step 2 "ECO" & "E":** Read points here.

**Step 5:** Chart the results. We provide a step-by-step manual with sample forms.

Step	Speed	Contrast	Relative Speed	Average Gradient
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				

## Specifications

System - microprocessor controlled.  
Film size reach - up to 7" (middle of 14" width film)  
Dimensions: 9.5" l x 3.5"h x 7.5"d (240 x 85 x 190mm)  
Weight: approx: ECO-2.2 lbs (1 kg) - E & I-3.3 lbs (1.5 kg)  
Power supply: 9.5 volt battery for short term use. Also provided with 120/220 power adapter

## Sensitometer

21-density steps - with increase of 0.15D per step  
Green color peak wavelength - 510nm  
Blue color peak wavelength - 460nm  
Exposure tolerance within each step - less than 0.01D  
Repeatability: better than +/-0.02 log e  
Time Stability: better than +/-0.02 log e  
Setting: separate adjustments for blue or green wavelength

## Autoreading Densitometer

Reads 21 steps in less than 5 seconds  
Automatic zero calibration  
No warm-up time required  
Information displayed on LCD  
Immediate evaluation of the following statistics:  
Speed Index & Contrast Index  
Base and Fog  
D step 21 (D-max)  
LE (Speed) & LK (Contrast)  
Relative Speed & Average Gradient  
Can display values for individual steps  
Can be set for English, Spanish, French and German  
RS232 Interface (for download into special software)  
Measure range: 1.00 - 4.50D  
Measure accuracy: +/- 0.01  
Stability: +/- 0.02D  
Independent light area - 3mm diameter  
Temperature operating range: 59 to 95F (15-35C)

## Point Densitometer

Length of measure arm - 8" (20cm)  
LCD display of results  
Calibration with enclosed calibration film  
Independent light area - 3mm diameter  
Measure range: 1.00 - 4.50D  
Measure accuracy: +/- 0.01  
Stability: +/- 0.01  
Temperature operating range: 59 to 95F (15-35C)

## Model Numbers

40450 - Densonorm 21-ECO  
40460 - Densonorm 21-E  
40470 - Densonorm 21i