

Elite Imaging Table



*Unobstructed 3-D
C-Arm Imaging*

500 lbs 227 kg Weight Capacity

Accessories



Wireless Hand Control



Foot Control for all Movements



Accessory Clamp



Patient Arm Board and Pad

0047-ELITE SERIES SPECIFICATIONS

		TM-5	TM-4	TM-3	TM-2	TL-2	TT-2	TE-1
(1) Motorized Elevation (Lift) Travel 300mm Lift Standard	Minimum	31" 79 cm	30" 76 cm	29" 73.5 cm	23" 58 cm	26" 66 cm	26" 66 cm	23" 58 cm
	Maximum	42.5" 108 cm	41.5" 106 cm	40.5" 103 cm	34.5" 88 cm	37.5" 95 cm	37.5" 95 cm	34.5" 88 cm
(1) Motorized Elevation (Lift) Travel 400mm Lift (Optional)	Minimum	NO	NO	NO	24.75" 630mm	28" 710mm	28" 710mm	24.5" 620mm
	Maximum	NO	NO	NO	44" 1020mm	43.5" 1150mm	43.5" 1150mm	40" 1020mm
(2) Motorized Longitudinal Table-top Travel		15¾" 400mm	15¾" 400mm	NO	15¾" 400mm	NO	NO	NO
(3) Motorized Lateral Table-top Tilt		± 15°	± 15°	± 15°	NO	± 15°	NO	NO
(4) Motorized Trendelenburg Tilt		± 15°	± 15°	± 15°	NO	NO	± 15°	NO
(5) Motorized Lateral Table-top Travel		± 4" 100mm	NO	NO	NO	NO	NO	NO
Patient Weight Capacity		500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg
Carbon Fiber Table-top 22" wide x 82" long, with Facial access cutout		YES	YES	YES	YES	YES	YES	YES
Carbon Fiber Rectanaglar Table-top 22" wide x 82" long, with Facial access cutout		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Un-obstructed Table-top Imaging Area		58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm
Radiolucent Sofcare Pad		YES	YES	YES	YES	YES	YES	YES
Pateint Restraints (3)		YES	YES	YES	YES	YES	YES	YES
OPTIONS								
Dual Patient Arm Boards with Pads		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Snap-on Arm Board with Pad and Mounting Clamp		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Wireless Hand Control with Receiver		NO	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Foot Pedel Control		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Rectangular Table Top		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL



North American Office: 23151 Alcalde Dr. Suite C-8 Laguna Hills, Calif. 92653
Worldwide Headquarters: Anavagen 1, Vaxjo SE-325 46, Sweden

High Quality Hospital Grade Mobile C-Arm Imaging Table

The “Elite” series of imaging tables are designed for high patient throughput utilizing the highest quality components and designed for ease of use by the medical staff with patient comfort in mind.

The table top has been designed to give maximum unobstructed Imaging allowing 360° rotation radiation free coverage around the longitudinal axis of the carbon fibre table-top which makes it ideal for 3D imaging studies.

The table-top pad is also designed to give the patient maximum comfort being manufactured utilizing our patient form fitting Visco/Memory “Sofcare” High Density foam.

150 mm (58”) Full 360° Unencumbered Rotational Imaging Area



Movement One

Motorized Elevating Lift

Lift Travel: 300 mm (11.75”)

Optional Travel: 400 mm (15.75”)

Movement Two:

Motorized Trendelenburg Tilt

+/- 15 Degree Tilt



Movement Three

Motorized Lateral Tilt

Lateral Tilt, +/- 15 Degree Tilt



Movement Four

Motorized Extend and Retract

Longitudinal Travel, (16") 40 mm



Movement Five

Motorized Lateral Drive

Lateral Travel, (+/- 4") +/-10 mm