



DIGITAL X-RAY

MOBILE SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/mobile



Wireless DR Panel Choices



 Cesium Scintillator


 Wireless Interface


 Cassette Sized


 14x17"

MOMENTUM WIRELESS FLATPANEL DR


2020imaging.net/momentum



 16-bit Images

 Lithium Ion Capacitor
(7-yr warranty)

 Featherlight 5.5 lbs


 Rapid Charge

VIZION+ DR MOBILE

WIRELESS FLATPANEL DR

2020imaging.net/vizion



 14-bit Images


 Lightweight 7.93 lbs



Hospital-Grade Software



ACQUISITION INTERFACE


- Under 6 second  Processing
- Preset X-ray Exams
- Unlimited Image Container
- Advanced Image Processing
- Auto-Shutter & Image Enhance
- Background Submission
- Adult & Pediatric Settings
- Fully Customizable
- Image Enhancement Tools
- User Friendly Interface

Multi-functon Studylist

- Auto DICOM Send Included
- Custom Search/Query
- Burn/Import CDs
- Customized User Settings

Advanced Viewer









- Full Image Viewer on the go
- Fast, User Friendly Interface
- Image Fine Tuning
- Customized Layout



Mobile SR-130 X-Ray



SR-130 FEATURES

-  High Frequency Inverter
-  Lower Radiation Dose*
-  Agile Lightweight Portability
-  Hospital Quality Imaging
-  Intuitive Interface
-  Superb Image Quality
-  Extended Tube Life
-  UL & CUL Listed
-  Built-in Cassette Storage

2020imaging.net/sr130

SR-130 Specifications

- 3.0kw Output - Regulated mA/kVp
- KVP Range: 40 - 100
- MA Range: 15 & 30
- MAS Range: 0.15 - 120





Durable
Superior
Efficient
Reliable
Affordable

2020imaging.net/mobile

