



WORLD-CLASS MEDICAL PRODUCTS NOW PACKAGED TO SHIP WORLDWIDE!

Model Number:


EVO-460 • Flat Pack Barrier

Key Features:

- Easy assembly: one box, one wrench and 4 bolts (all included).
- Lower freight cost.
- Lightweight, super-high-strength aluminum frame.
- Hospital-grade casters for easy mobility, 2 with brakes.
- .50mm LE top to bottom.

Details:

Dimensions	72" H x 34½" W x 25" D
Window Size	24" H x 30" W
Protection	Lead free acrylic window .50mm LE
Frame	Aluminum
Panel	.50mm LE in White, Black, or Grey
Wheels	4 swivel hospital grade casters (2 with brakes)
Floor Curtain	Removable .50mm LE curtain included



Packaging:



Overall Dimensions	11½" H x 40½" W x 50½" D
Total Weight	135 lbs.

* 5 piece set • Easy to assemble with everything you need included!

