## SAFETY DATA SHEET

## SECTION 1 - General Information PRODUCT NAME Med-Dent TM Developer Premixed Health 1 С Photographic Solution HMIS SCALE: Flamability 0 Reactivity 0 PPE PRODUCT USE White Mountain Imaging CATALOG #: SUPPLIER PO Box 216 Salisbury, NH 03268 NS-204, NS-324 ADDRESS 603-648-2124 EMERGENCY TEL# CHEMTREC 800-424-9300 PHONE NO SECTION 2 - Hazards Identification Precautionary Statement. Category <u>GHS</u> WARNING OSHA/GHS Oral: May be harmful if swallowed 5 n/c 5 Skin: May cause mild skin irritation or reaction. n/c Eye: May cause eye irritation 2B 2B EFFECTS OF CRONIC EXPOSURE TO PRODUCT None known. (n/c: not within OSHA classification limits.) SECTION 3 - Product & Hazardous Ingredients Information **HAZARDOUS INGREDIENTS CAS NUMBER** W/W % Permissible Exposure Limit LD50 7732-18-5 Water 75-85 Not listed 42.6gm/kg Sodium metabisulfite NOIOSH TWA 5 mg/m3 2000 mg/kg 7681-57-4 < 4 375 mg/kg 123-31-9 OSHA TWA 2 mg/m3 Hvdroquinone < 2 Potassium hydroxide 1310-58-3 < 3 NIOSH © 2mg/m3 365 mg/kg **SECTION 4 - First Aid Measures** Wash exposed skin thoroughly with soap and water. If irritation persists seek medical advice. EYE CONTACT Flush eyes with water for 15 minutes. Seek medical advice if symptoms persist. DO NOT induce vomiting. Seek medical attention immediately. INHALATION If inhaled, remove to fresh air. SECTION 5 - Fire Fighting Measures YES NO IF YES UNDER WHAT CONDITIONS Flamability MEANS OF EXTINGUISHING As needed to extinguish adjacent fire source. FLASH POINT (C/F) AND METHOD UPPER LIMIT ( % BY VOLUME) LOWER LIMIT ( % BY VOLUME) AUTOIGNITION TEMPERATURE (C/F) Not appropriate Not appropriate Not appropriate Not appropriate HAZARDOUS COMBUSTION PRODUCTS EXPLOSION DATA SENSITIVITY TO IMPACT SENSITIVITY TO STATIC DISCHARGE "See Hazardous Decomposition Products" Not susceptable Not susceptable **SECTION 6 - Accidental Release** Flush to sewer with plenty of water, if permitted. If not, soak up with dry absorbent material. Consult with governmental regulatory agencies for appropriate disposal of this material. Wear appropriate protective equipment. SECTION 7 - Handling and Storage

Read and follow the label, safety data sheet, and instructions before using. Avoid any contact with skin and eyes. Keep container tightly closed. Avoid incompatible substances. Wash thoroughly after handling this product.

PREPARED BY

		SECTION 8	3 - Exposure Controls			
VENTILLATION Good general vent skin Latex or nitrile			re: Splash goggles	RESPIRATORY	Not normally required	
OTHER	HER Eye wash station					
		SECTION 9 - Phy	sical & Chemical Properti	es		
PHYSICAL STATE Liquid		ODOR & APPEARANCE Light red, odorless solution.			odor threshold  Not available	
VAPOR PRESSURE (mm-Hg) <17 mm Hg		vapor density (Air=1) 0.6 mm Hg	EVAPORATION RATE (Water=1)  Not Available	BOILING POINT > 212 F	FREEZING POINT < 32 F	
рН 10.20		SPECIFIC GRAVITY 1.08	SOLUBILITY IN WATER  Complete	COEFF WATER/OIL N/A	%VOLATOILE (BY VOLUME) 80-90% as water	
		SECTION 10 - S	Stability & Reactivity Data			
CHEMICAL STABILIT	YES X	IF NO, UNDER	WHAT CONDITIONS			
INCOMPATABLE WIT	TH OTHER SUBSTANCES		ONES? ing agent and strong acids.			
REACTIVITY AND UN	IDER WHAT CONDITIONS	Stable under	ambient pressure and temp	erature		
HAZARDOUS DECO	MPOSITION PRODUCTS	None	ambioni procedio and tomp	orataro		
		SECTION 11 - 7	Foxicological In formation			
EXPOSURE LIMITS		IRRITANCY OF PRODUCT	SENSITIZATION TO PRODUCT		CARCINOGENICITY	
None known		Skin: Irritant	None known		Not listed by IARC, OSHA, or ACGIH	
TERATOGENICITY		REPRODUCTIVE TOXICITY MUTAGENIC			SYNERGISTIC PRODUCTS	
None Known		Not Known	None known		None known	
		SECTION 12	- Ecological Information			
No Data Avai	ilable					
		SECTION 13 -	Disposal Considerations			
Dispose of th	is product in acc	cordance with federal, state a	and/or local laws.			
		SECTION 14 - T	ransportation Information	1		
Not regulated	d for transportation	on.				
		SECTION 15	- Regulatory Information			
SARA Sectio SARA Sectio	n 304, CERCLA n 302,304,313 h	RQ) for Hydroquinone 5000 l . Hazardous components - (Nazardous components (Non Hazardous components (Non er Canadian WHMIS regulati	None) e) ons.	0 lb		
		SECTION <sup>2</sup>	16 - Other Information			
		his SDS is furnished without vered and must make an inde				

White Mountain Imaging PHONE NO 603-648-2124 Date: 5/20/2015

Cat#: NS-2014, NS-324

assure proper use and disposal of this material and the health and safety of employees and customers.