SAFETY DATA SHEET

SECTION 1 - General Information PRODUCT NAME Med-Dent ™ Fixer Premixed Photographic Solution C HMIS SCALE: Health 1 Flamability 0 Reactivity 0 PPE PRODUCT USE White Mountain Imaging CATALOG #: SUPPLIER PO Box 216 Salisbury, NH 03268 NS-204, NS-324 ADDRESS 603-648-2124 CHEMTREC 800-424-9300 **EMERGENCY TEL#** PHONE NO **SECTION 2 - Hazards Identification** Category WARNING Precautionary Statement. OSHA/GHS GHS 5 Skin: May cause mild skin irritation and or reaction. n/c Eye: Can cause irritation of eyes and mucous membranes. 2B 2B n/c 5 Oral: May be harmful if swallowed EFFECTS OF CRONIC EXPOSURE TO PRODUCT None known. (n/c: not within OSHA classification limits.) SECTION 3 - Product & Hazardous Ingredients Information **HAZARDOUS INGREDIENTS** Permissible Exposure Limit **CAS NUMBER** W/W % LD50 Water 7732-18-5 90-95 Not listed 42.6 gm/kg 2890 mg/kg Ammonium Thiosulfate 7783-18-8 8-12 Not listed Sodium Sulfite Not listed 7757-83-7 2825 mg/kg <2 Sodium Acetate 127-09-03 Not listed 3530 mg/kg **SECTION 4 - First Aid Measures** SKIN CONTACT Wash exposed skin thoroughly with soap and water. If irritation persists seek medical advice. Flush eyes with water for 15 minutes. Seek medical advice if symptoms persist. DO NOT induce vomiting. Seek medical attention immediately. INHAI ATION If inhaled, remove to fresh air. **SECTION 5 - Fire Fighting Measures** Flamability YES NO IF YES UNDER WHAT CONDITIONS MEANS OF EXTINGUISHING As needed to extinguish adjacent fire source. FLASH POINT (C/F) AND METHOD UPPER LIMIT (% BY VOLUME) LOWER LIMIT (% BY VOLUME) AUTOIGNITION TEMPERATURE (C/F) Not appropriate Not appropriate Not appropriate Not appropriate HAZARDOUS COMBUSTION PRODUCTS EXPLOSION DATA SENSITIVITY TO IMPACT SENSITIVITY TO STATIC DISCHARGE "See Hazardous Decomposition Products" Not susceptable Not susceptable **SECTION 6 - Accidental Release** Flush to sewer with plenty of water, if permitted. If not, soak up with dry absorbent material. Consult with governmental regulatory agencies for appropriate disposal of this material. Wear appropriate protective equipment.

SECTION 7 - Handling and Storage

Read and follow the label, safety data sheet, and instructions before using. Avoid any contact with skin and eyes. Keep container tightly closed. Avoid incompatible substances. Wash thoroughly after handling this product.

ABX_{TM} Fixer Premixed

		SECTION	8 - Exposure Controls		
VENTILLATION SKIN OTHER	Good general ventillation Latex or nitrile Eye wash station		EYE: Splash goggles RESPIRATORY		Not normally required
		SECTION 9 - Ph	ysical & Chemical Properti	es	
PHYSICAL STATE ODOR & APPEARANCE Light Blue liquid with		ODOR & APPEARANCE Light Blue liquid with slight	slight vinegar odor.		odor threshold Not available
VAPOR PRESSURE (mm-Hg) <17 mm Hg		VAPOR DENSITY (Air=1) 0.6 mm Hg	EVAPORATION RATE (Water=1) Not Available	BOILING POINT > 212 F	FREEZING POINT < 32 F
рН 4. 4	40	SPECIFIC GRAVITY 1.085	SOLUBILITY IN WATER Complete	COEFF WATER/OIL N/A	%VOLATOILE (BY VOLUME) 80-90% as water
		SECTION 10 -	- Stability & Reactivity Data		
CHEMICAL STABILITY	YES X	IF NO, UNDER	R WHAT CONDITIONS		
INCOMPATABLE WITH	HOTHER SUBSTANCE		CH ONES? ucing agents & Strong bases.		
REACTIVITY AND UNI	DER WHAT CONDITION	Stable unde	er ambient pressure and temp	perature	
HAZARDOUS DECOM	MPOSITION PRODUCTS		kide & ammonia		
		SECTION 11 -	· Toxicological In formation	<u> </u>	
		IRRITANCY OF PRODUCT	SENSITIZATION TO PRODUCT		CARCINOGENICITY
None known		Skin: Irritant	None known		Not listed by IARC, OSHA, or ACGIH
TERATOGENICITY		REPRODUCTIVE TOXICITY	MUTAGENICITY	MUTAGENICITY	
None Known		Not Known None kno		n	None known
		SECTION 12	2 - Ecological Information		
No Data Avail	lable				
		SECTION 13	- Disposal Considerations		
Dispose of thi	is product in ac	cordance with federal, state	and/or local laws.		
		SECTION 14 -	Transportation Information	า	
Not regulated	for transportat	ion.			
		SECTION 15	5 - Regulatory Information		
SARA Section SARA Section	n 304, CERCLA n 302,304,313	RQ) for Ammonium Thiosulf A Hazardous components - A Hazardous components (No ler Canadian WHMIS regula	Ammonium Thiosulfate one) tion		
		SECTION	l 16 - Other Information		

The information contained in this SDS is furnished without warranty of any kind. The user should consider this data a supplement to other data gathered and must make an independent determination of suitability from this and other sources to assure proper use and disposal of this material and the health and safety of employees and customers.

> PHONE NO 603-648-2124 Date: 5/20/2015 NS-204, NS-324 Page 2/2