Porta 100HF for Veterinary TELEVERE

"The Green Machine - Dual 100HF" High Frequency Portable X-ray Unit





100 kVp, 30MA Dual Distance Lasers 26" SID Weighs Only 19 lbs (Lightest 100 kVp on the market) Six (6) Year Parts Warranty Perfect X-ray Machine for All Veterinarians: Small Animal, **Equine and Mixed Animal**









Porta 100HF for Veterinary

"The Green Machine"

High Frequency Portable X-ray Unit



Specifications

"The Green Machine™" Porta 100HF is the newest technology to hit the Veterinary market. It is the only portable x-ray machine on the market boasting a (6) six year parts warranty. Manufactured in Japan with the highest quality parts available. This powerful and light x-ray machine weighs only 19 lbs.

"The Green Machine™" is versatile and perfect for small animal, equine, and mixed animal practices.

System Features

- High Frequency High Output 2 kW, 100 kVp, 30 mA
- Rotary Switches for easy and quick selection of kVp and mAs
- Programmable APR stores up to (8) presets
- Tape measure
- DR Flat Panel compatibility with modular sync connector
- Automatic mAs control for line voltage compensation (helpful with long extension cords)
- Dual Distance lasers preset to 26" providing accurate FFD for Equine Usage
- Aluminum Carrying Case for easy transport
- (6) Year Standard Parts Warranty

Product Specifications

- kVp setting: 40-100 kV (in 2 kVp step)
 Rotary switch for digital display with 7 segment LED
- mAs selection: 0.3-20 mAs (32 settings)
- X-ray tube:
 - a) Focal spot: 1.2 x 1.2 mm
 - b) Anode heat storage capacity: 20 kHU
 - c) Target angle: 16 degrees
- Power requirement: 100-120 V or 200-240 V
- Weight: 19.4 lbs

Upgrade Options

- Maximal View X-ray Table
- Digital DR Systems Equine
- Digital DR Systems Small Animal
- Digital CR Systems

