

6-Drawer X-Tall Emergency Cart (UXRLA-333669-RED)

 waterloohealthcare.com/product/aluminum-emergency-carts







Waterloo made its name in emergency carts back in the 1950s, and is still the most recognized name in emergency carts today. Our elite aluminum carts are lighter and more advanced than others allowing for a fast and efficient code response.

KEY FEATURES:

- 30% lighter than steel, for maximum speed and efficiency
- Dual pull out shelves, doubling your work surface
- Adjustable metal dividers included in each drawer
- Transport lock keeps drawers secured during rapid transit

SKU: UXRLA-333669-RED Category: Emergency Carts Tag: Emergency Cart Aluminum

SPECIFICATIONS

Frame Type: ALUMINUM
STANDARD
Frame Size: X-TALL
Dimensions: 29"w x 24.5"d x 45"h




Drawer
Dimensions: 22"w
x 16.5"d



Drawer Height: 3-3-3-6-6-9

3"		3 drawers
6"		2 drawer
9"		1 drawer
12"		

Frame Color: RED



Drawer Color:
RED



Lock System: LEVER LOCK



Top: REMOVABLE PLASTIC



Bumpers:
STABILIZER
BASE



Casters: 5" STANDARD

