

## 4-Drawer Short Emergency Cart (USBLA-3369-ELB)

 [waterloohealthcare.com/product/4-drawer-short-usbla-3369-elb](https://waterloohealthcare.com/product/4-drawer-short-usbla-3369-elb)







Waterloo made its name in emergency carts back in the 1950s, and is still the most recognized name in emergency carts today. Our elite aluminum carts are lighter and more advanced than others allowing for a fast and efficient code response.

#### KEY FEATURES:

- 30% lighter than steel, for maximum speed and efficiency
- Dual pull out shelves, doubling your work surface
- Adjustable metal dividers included in each drawer
- Transport lock keeps drawers secured during rapid transit

SKU: USBLA-3369-ELB Category: Emergency Carts Tag: Emergency Cart Aluminum


## SPECIFICATIONS

---

Frame Type: ALUMINUM  
STANDARD  
Frame Size: SHORT  
Dimensions: 29"w x 24.5"d x 36"h



Drawer  
Dimensions: 22"w  
x 16.5"d



Drawer Height: 3-3-6-9  
2  
drawers



1  
drawer  
  
1  
drawer


Frame Color: RED




Drawer Color:  
RED




Lock System: LEVER LOCK



Top: REMOVABLE PLASTIC



Bumpers:  
STABILIZER  
BASE



Casters: 5" STANDARD

