



Pedia-Poser

The *Pedia-Poser* provides fast, gentle, and secure immobilization of infants through 4-year-olds for X-ray imaging. Its uniquely designed saddle seat, with detachable saddle horn, is form-fitting and very safe. The bench seat for older children includes a detachable saddle horn that aids in positioning.

- Immobilization is done gently, but securely
- Eliminates virtually all retakes, reducing radiation exposure, as well as saving time and money
- The patient is secured only once – the chair rotates, positioning for different views
- Keeps patient as close to the imaging detector as possible to ensure high-quality images
- Adjustable back and head-tilt support fits patients of all sizes
- Critical chair components are X-ray lucent for excellent X-ray contrast
- Easy to keep clean, the *Pedia-Poser* is made from germ-resistant, non-porous polyethylene, and stainless steel
- High-quality, lockable casters provide mobility and stability
- Simple, reliable, and virtually indestructible

**MADE IN THE
USA**

Clear Image Devices ^{LLC}
5225 68th St., Suite 8 Caledonia, MI 49316
Phone: +1 734-474-6537
Fax: +1 734.661.0220



PEDIA-POSER			
Item #	Description		
24301	For infants through 4-year-olds One replacement strap pack included CID Practice Positioning Bear included		
Image Type	AP	Lateral	Oblique
C-Spine	X	X	X
Chest	X	X	
Abdominal	X		
Magnified Airway	X	X	
Dimensions:			
Infant Seat Width:	8.4 in / 21.3 cm		
Toddler Seat Width:	9.5 in / 24.1 cm		
Wheel Base Width:	27 in / 68.6 cm		
Seat Height:	25.25 in / 64.1 cm		
Total Chair Height:	61.75 in / 156.8 cm		
50123	Additional replacement strap pack		

Free Standard Shipping - All Products - 48 States

Clear Image Devices has become the #1 brand in imaging accessories through innovative product development in cooperation with the medical imaging community. We believe that value begins with quality; quality products and customer service. We are advancing medical imaging every day by bridging the gap between technology and people.

www.ClearImageDevices.com

copyright ©2016 Clear Image Devices LLC
All rights reserved.