



RADIATION GOGGLES & FACESHIELDS

MODEL RG-FM-FULL



RG-FM-FULL

DESCRIPTION

Phillips Model FM-FULL is a radiation protective mask that eliminates blind spots caused by standard glasses and offers a wide, unobstructed viewing area. The protective visor has a scratch-free coating and offers .10mm lead equivalency. The mask includes a ratcheting head strap that fits all sizes and is especially useful for people who wear glasses or need to protect their entire face area.

FEATURES

- *0.10mm Lead Equivalency Radiation Protective Mask*
- *Ratcheting Head Strap to Fit All Head Sizes*
- *Scratch Coating on Lead Acrylic*

SPECS

Frame Size: ADJ
Weight: 460g
CE Pending
Non-Rx Frame