

Article Information Sheet

Producer: FUJIFILM Corporation
MIDTOWN WEST, 7-3, AKASAKA 9-CHOME, MINATO-KU, TOKYO 107-0052, JAPAN
Distribution: FUJIFILM Medical Systems U.S.A., Inc.; 419 West Avenue, Stamford, CT 06902;
Tel.: (203) 3242000, Fax: (203) 3530926;
Issued: 4/7/1997
Modified: 4/1/1998

1. Product

Product Name: FUJI MEDICAL IMAGING PLATE

2. Composition

All products in this category consist of a coated polyester support.
The photosensitive layer is fixed mainly by resin.
The main ingredients are:

<u>Main Ingredients:</u>	Weight-%:
polyester	30-50
barium fluorhalide	40-60
resin	8-12
additives	1-5

3. Health & Safety Hazards

If used as intended, this product does not have any special health or safety concerns.

4. Fire Hazards and Fire Fighting Measures

This product is capable of burning, but will not ignite or flash under normal conditions.

<u>Product / Main Ingredient(s):</u>	Ignition Temp.	Flash Point	Density
product	n/av	n/av	n/av
polyester	n/av	n/av	n/av

Combustion Gases:

Carbon dioxide, carbon monoxide, hydrogen halides and small amounts of nitrogen oxide can be generated depending upon the burning condition.

Fire Fighting Measures:

Water, carbon dioxide, extinguishing powder or foam can be used as an extinguishing media.

While fire fighting, wear protective equipment such as self contained breathing apparatus, depending upon the fire situation.

5. Disposal Considerations

5.1 General Disposal Considerations

It is recommended to separate the waste plates from other waste streams and to provide them to a legitimate recycler.

5.2 Special Disposal Information for the US

6. Transport

This product is not subject to dangerous material classification and labelling pertaining to transportation.

7. Storage

Specific hazards regarding storage are not known. For information on storage conditions please refer to the technical data sheets.

8. Further Information

Further information is not available at present.

The content of this information sheet has been drawn up in accordance with current knowledge at the time when the sheet was issued. This Article Information Sheet provides information related to the environment and to safety and does not make any warranty about the properties or quality of the product. Furthermore, the precautions presented herein are postulated upon normal handling (the type of use originally intended by the manufacturer) of the product.