

SPECIFICATIONS

Immobilizer is approximately 20" x 20" x 35" (50 cm x 55 cm x 88 cm) in height. Stainless steel table with 3" (7.62 cm) casters, 2 are locking. Table top and shelf are made of extremely durable ABS composite material in child-pleasing primary colors which requires no maintenance and is virtually indestructible. All other moving parts are heavy plexiglass.

GUARANTEE LIMITED

Modern Way Immobilizers, Inc. guarantees the Pigg-O-Stat Immobilizer to the extent that any part or parts returned to us postpaid within one year from date of purchase and found to be defective in material or workmanship by our inspection will be replaced free of charge, fob our factory.

Even though the plexiglass parts will withstand many times the stress and strain required to perform their purpose under proper use, they are not guaranteed against breakage since the manufacturers of plexiglass will not make such guarantee to us.

Under no circumstance will Modern Way Immobilizers, Inc. be responsible for misuse of this equipment or for any guarantee other than the above context made and signed by the president of Modern Way Immobilizers, Inc.

Modern Way Immobilizers, Inc.

Jimmy C. Pigg

Jimmy C. Pigg, President

Jalmer Pigg, Sr.
1917 - 1995
Founder



MODERN WAY IMMOBILIZERS, INC.

Producing Quality Products Since 1960

P.O. BOX 660 • 104 JOHNSON STREET • CLIFTON, TN 38425

PHONE 931-676-5274 • FAX 931-676-3860

TOLL FREE 866-MWI-PIGG or 866-694-7444

www.piggostat.com



U.S. PATENT NO. 5,600,702

PIGG-O-STAT[®]

IMMOBILIZER and POSITIONER

EASILY ADJUSTS

ROTATES | FITS FORM | PRECISE

