

# HORIZON RETRACTO™ RAIL HYDRAULIC STRETCHER



DESIGNED FOR PROCEDURE, RECOVERY & TRANSPORT

## 462 SERIES RETRACTO™ RAIL HORIZON HYDRAULIC STRETCHER

The Retracto™ Rail Hydraulic Stretcher is designed to provide longlasting and consistent performance. These easy-to-clean units support patients during transport in any hospital or surgical environment.

### VERSATILITY:

- Streamline patient handling allowing this unit to be used from preinduction to transport to recovery to help improve case turnaround
- Easy to use hydraulic foot pedals to adjust height and Trendelenburg on either side of the stretcher
- Intermediate side rail position offers maximum patient access without compromising patient or staff safety
- Dual accessory wells at all four corners for versatility in adapting accessories
- Full perimeter heavy duty non-marring bumper system
- Seamless thermoformed base cover with storage
- Built-in oxygen tank holder



### Specifications:

	inch	cm†
Overall Length:	83	210
Overall Height: <i>range - minimum</i>	22.25	56
Overall Height: <i>range - maximum</i>	34.75	88
side rail Length:	60	152
side rail Height:	15	38
Tente® Casters: <i>easy rolling</i>	8	20
	degree of movement †	
Cable Free Fowler: <i>pneumatic-assist</i>	0 - 90	
Trendelenburg: <i>standard</i>	0 - 16	
Trendelenburg: <i>reverse</i>	0 - 16	
	lb	kg†
Shipping Weight:	330	150

### MOBILITY:

- Three modes of operation: steer, neutral and brake
- Brake and steer levers are accessible at all four corners of the stretcher base
- 8 inch Tente® Casters

### Features:

Rail Type:	Retracto™
IV Wells:	4 - 1/corner
Accessory Wells:	4 - 1/corner
Brake/Steer Pedals:	4 - 1/corner
Hydraulic Controls:	Foot Pedal - 1/side

## 426-HMC-ST HYDRAULIC MODEL

*shown with back tilted up*

 462 Product Info



	462-EMC-ST		462-FMC-ST		462-GMC-ST		462-HMC-ST	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Overall Width:	30	76	34.5	88	30	76	34.5	88
patient surface Width:	25	64	30	76	26	64	30	76
patient surface Length:	76	193	76	193	76	193	76	193
	feature		feature		feature		feature	
Knee Gatch:	No		No		Yes		Yes	
	lb	kg†	lb	kg†	lb	kg†	lb	kg†
Weight Capacity: EVENLY DISTRIBUTED	500	227	800	363	500	227	800	363

■ Standard ■ Optional — Does Not Apply  
† same specs for both AC110 volt and the DC220 volt versions  
‡ standard rounding method used to show metric values as whole numbers

Product Hyperlinks Enabled