

# Product Specification Sheet

## PRODUCT TITLE:

Specification of Product 1:	<b>Bind-It™ Radioactive Decontaminant Concentrate</b>
Specification of Product 2:	<b>Bind-It™ Radioactive Decontaminant Ready-to-Use Spray</b>
Specification of Product 3:	<b>Bind-It™ Radioactive Decontamination Hand and Body Soap</b>
Specification of Product 41:	<b>Bind-It™ Radioactive Decontaminant Surface Wipes</b>

Bind-It™ Decontamination products offer highly effective, non-corrosive, radioactive removal products that are pleasantly scented and safe for routine use. Bind-It is designed to be effective on most Isotopes and is especially effective on all radio-iodines.

Bind-It has been independently tested and is used by hundreds of nuclear pharmacies, nuclear medicine departments and veterinary hospitals. Bind-It is effective on I123, I125, I131, I129, Tc99m, TI201, In111, F18, Ga67, Mo99, Ac225, Cs137, Pb212 and more.

Bind-It™ Decontaminate is shelf stable, the stability is unaffected exposure to air. Bind-It™ Decontaminate has a guaranteed shelf life of at least 24 months, open or unopened.

## Key Features of Bind-it Decontamination Fluid

- 🧼 Effective cleaner and decontaminant in one.
- 🧼 Highly effective on: Tc, Mo, I, FI, TI, Ga, In, Cr, Sm - all NM Isotopes.
- 🧼 Pleasantly scented - no harsh odors!
- 🧼 Safe for use on syringe shields; will NOT cloud lead glass.
- 🧼 Non-corrosive, non-toxic, non-flammable with a pH of 6.5.
- 🧼 Safe on delicate detectors, probes & survey meters, all surfaces.
- 🧼 Cleans washrooms, bench tops, fume hoods, floors—most surfaces.
- 🧼 Extremely effective on I131; clean isolation rooms & patient “accidents”.

## Key Features of Bind-it Hand and Body Soap

- ☑ Gentle formula will not dry out your skin
- ☑ Pleasant almond scent
- ☑ Replaces harsh general purpose lab soaps
- ☑ Binds and removes radioactive skin contamination
- ☑ Economical

