

THE MOST PROTECTIVE APRON IN THE WORLD

Revolution Aprons have been designed to be the most protective x-ray protection aprons in the industry. Current style aprons only offer your desired protection levels at the overlap; however Revolution aprons protect you at your desired lead equivalency throughout the entire apron.

Example: A user desires a .5 mm Pb Eq apron, current styles offer you two .25 panels to achieve .5 mm Pb Eq. at the overlap. This leaves the shoulders and side area only protected at .25 mm Pb.



SHOULDER

Each panel of a Revolution Vest & Skirt is .5 mm Pb. Compared to traditional styles of aprons that use 2 panels of .25 mm to get to .5 mm Pb

PARTIAL OVERLAP

Revolution Aprons offer you the desired protection levels throughout, without relying on a full overlap.

THE MOST COMFORTABLE APRON IN THE WORLD

The Revolution Vest and Skirt is the most advanced x-ray protection apron on the market today. We designed it to blend the highest levels of protection with comfort and style and have achieved that goal. Revolution Vest and Skirts are truly the most comfortable aprons in the world.



BACK FLAP DESIGN

Revolution Vests feature our all new "Back Flap" design, providing extra protection when bending at the waist as well as balancing out the weight of the apron front —to-back for increased comfort.

BUILT-IN BLACK BELT

Revolution Skirts have been designed with a built-in orthopedic elastic Black Belt. This added support relieves weight on the lower back and provides a superior fit for the skirt. The skirt waist is now adjustable depending on how much support you would like on your lower back. The built in Black Belt also stops the skirt from slipping down the waist which provides a much more comfortable and secure fit while completely eliminating the need for aids like suspenders or hot and uncomfortable "non-slip" materials.



View Online
InfabCorp.com

