

Nothing protects like

Wolf

Mobile Protection

The Clear Lead Mobile Barrier provides a full-body shield for one or two doctors or radiology technologists. The barrier can be easily repositioned or moved from room to room. The optically clear, lead plastic viewing panel rests on an opaque lower panel, and provides exceptional visibility with 0.5mm or 1.0mm lead-equivalent protection. The shatter resistance and durability of Clear Lead ensures long life under heavy use.



CLEAR LEAD BARRIER SPECIFICATIONS

Overall Size: 75" H x 31" W x 25" D	PB
56603 Lead Plastic Panel: 24H x 30W	0.5MM
56606 Lead Plastic Panel: 48H x 30W	0.5 MM
56604 Lead Plastic Panel: 54.5H x 30W	0.5MM
Overall Size: 75"H x 49"W x 25"D	PB
56600 Lead Plastic Panel: 30H x 48 W	0.5MM
56605 Lead Plastic Panel: 36H x 48W	1.0MM

Wolf
X-RAY CORPORATION

100 West Industry Court
Deer Park New York 11729

www.wolfxray.com
800-356-9729

