

DELWORKS

Workstation

High resolution viewing, intuitive interface

Features

Large 23" Widescreen, 2MP Multi-touch Color LCD Monitor

Powerful 3.0GHz Processor

16GB RAM, 1TB HD



Specifications

Computer Hardware Technical Specifications

Operating System	Windows 7 Pro 64-bit
Processor	3.0GHz, Intel Core i5
Storage/Memory	16GB RAM, 1TB HD

Monitor Technical Specifications

Cabinet Color	Black, Gray
Panel Type	Color TFT LCD Panel (IPS)
Panel Size	58 cm / 23" (584 mm diagonal)
Display Size	509.1 mm x 286.4 mm
Pixel Pitch	0.270 x 0.270 mm
Native Resolution	1920 x 1080
Color Support	16.77 million from a palette of 1.06 billion colors
Viewing Angles	178°, 178°
Brightness	260 cd/m ²
Contrast Ratio	1000:1
Respose Time	16 ms (On/Off)
Dimensions	21.9" x 14.4" x 8.1" (W x D x H)
Net Weight	14.55 lbs-with stand

08.2015 #9000-DW-Workstation



www.delmedical.com



New York
28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Chicago
241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

DEL MEDICAL
A **UMG** COMPANY