



FEMALE



MALE



Female Frames

	X-Small	Small	Medium	Large	X-Large	XX-Large
Bust	30-32" / 76-81 cm	33-36" / 84-94 cm	37-39" / 94-99 cm	40-42" / 102-107 cm	43-45" / 109-114 cm	46-48" / 117-122 cm
Waist	25-27" / 64-69 cm	28-30" / 71-76 cm	31-33" / 78-84 cm	34-36" / 86-91 cm	37-39" / 94-99 cm	40-42" / 102-107 cm

Male Frames

	Small	Medium	Large	X-Large	XX-Large
Chest	34-37" / 86-94 cm	38-41" / 97-104 cm	42-45" / 107-114 cm	46-49" / 117-125 cm	50-52" / 127-132 cm
Waist	28-31" / 71-79 cm	32-35" / 81-89 cm	36-39" / 91-99 cm	40-43" / 102-109 cm	44-46" / 112-117 cm

Female Standard Apron Lengths

	Frontal Apron	Full Apron	Classic & 703 Vest	Classic Skirt	Revolution Vest Front	Revolution Vest Back	Revolution Skirt
X-Small	38" / 97 cm	38" / 97 cm	20" / 51 cm	18" / 46 cm	22" / 56 cm	26" / 66 cm	22" / 56 cm
Small	38" / 97 cm	38" / 97 cm	22" / 56 cm	20" / 51 cm	22" / 56 cm	26" / 66 cm	22" / 56 cm
Medium	38" / 97 cm	38" / 97 cm	24" / 61 cm	22" / 56 cm	22" / 56 cm	26" / 66 cm	22" / 56 cm
Large	38" / 97 cm	38" / 97 cm	25" / 64 cm	23" / 58 cm	22" / 56 cm	26" / 66 cm	22" / 56 cm
X-Large	38" / 97 cm	38" / 97 cm	26" / 66 cm	24" / 61 cm	22" / 56 cm	26" / 66 cm	22" / 56 cm
XX-Large	38" / 97 cm	38" / 97 cm	27" / 69 cm	25" / 64 cm	24" / 61 cm	28" / 66 cm	22" / 56 cm

Male Standard Apron Lengths

	Frontal Apron	Full Apron	Classic & 703 Vest	Classic Skirt	Revolution Vest Front	Revolution Vest Back	Revolution Skirt
Small	38" / 97 cm	38" / 97 cm	22" / 56 cm	20" / 51 cm	24" / 61 cm	28" / 71 cm	22" / 56 cm
Medium	38" / 97 cm	38" / 97 cm	24" / 61 cm	22" / 56 cm	24" / 61 cm	28" / 71 cm	22" / 56 cm
Large	38" / 97 cm	38" / 97 cm	25" / 64 cm	23" / 58 cm	24" / 61 cm	28" / 71 cm	22" / 56 cm
X-Large	42" / 107 cm	38" / 97 cm	26" / 66 cm	24" / 61 cm	24" / 61 cm	28" / 71 cm	22" / 56 cm
XX-Large	44" / 112 cm	38" / 97 cm	28" / 71 cm	25" / 64 cm	26" / 66 cm	30" / 71 cm	22" / 56 cm