

Radiation Resistant Glove Comparison Chart

Model	ASG	Attenuator X	IBG	<u>xguard</u> RR1 69874	<u>xguard</u> RR1 67159	<u>XGUARD</u> RR2 62134	<u>XGUARD</u> RR2 67258	<u>Secure</u> Touch XR1	<u>Secure</u> Touch XR2	<u>Revolution</u> Light	<u>Revolution</u> Medium	<u>Revolution</u> Heavy
Thickness	.23 mm	.35 mm	.18 mm	.22 mm	.22 mm	.30 mm	.30 mm	.20 mm	.35 mm	.20 mm	.32 mm	.42 mm
Length	11"	11"	11"	12"	12"	12"	12"	11"	11"	11.5"	11.5"	11.5"
Attenuation @ 60 KVP	58%	66%	56.00%	45%	45%	55%	55%	46%	64%	40%	68%	75%
Attenuation @ 80 KVP	47%	55%	47%	35%	35%	45%	45%	36%	54%	35%	61%	69%
Attenuation @ 100 KVP	40%	43%	41%	26%	26%	35%	35%	30%	48%	30%	54%	63%
Lead Free	yes	yes	yes	no	yes	no	yes	yes	yes	no	no	no
Latex Free	yes	no	yes	no	yes	no	yes	no	no	yes	yes	yes
Powder Free	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes
Price Per Pair	\$40.00	\$41.00	\$45.00	\$39.20	\$53.20	\$43.40	\$57.40	\$36.40	\$36.40	\$38.00	\$42.00	\$46.00
Available Sizes	6.5 - 9	6.5 - 9	6.5 - 9	6 - 9	6 - 9	6 - 9	6 - 9	6.5 - 9	6.5 - 9	6.5 - 9	6.5 - 9	6.5 - 9