

FABRIC

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables and passes testing for Cytotoxicity, Skin Irritation and Sensitization as well as CA TB117 Fire Retardance testing. For disinfecting purposes you may use Protex, MadaCide, MadaWipe, Accell TB, Virox® (Canada), PDI Super Sani-Cloth Germicidal wipe or a diluted 1/100 bleach solution.

TERRATOUGH™



WOOD FINISHES

Our water based finishes are highly scratch and water resistant. They can be wiped down with disinfectants as needed. Stains allow the natural grain of the wood to show through and create depth. Lacquers are more opaque and have a more solid appearance. **Available on select models.**

Standard
Natural Finish



Optional Stains



Optional Lacquers



FABRIC

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables and passes testing for Cytotoxicity, Skin Irritation and Sensitization as well as CA TB117 Fire Retardance testing. For disinfecting purposes you may use Protex, MadaCide, MadaWipe, Accell TB, Virox® (Canada), PDI Super Sani-Cloth Germicidal wipe or a diluted 1/100 bleach solution.

TERRATOUGH™



THERMOFORM

Choose from 4 high quality acrylic colors to match your decor. Available on Procedure Chairs only.



POWDER COATING

Powder coating colors are product specific with no options available. These are provided for reference only.

